

GUEDEL AIRWAY

oropharyngeal airway with ET Tube fixation possibility

- rail and lateral clamp design allow direct fixation of the ET Tube to the Guedel Airway
- hard biteblock but soft and atraumatic tube
- same lumen from proximal to distal end allows the use of largest possible suction catheters
- one side open with smooth inner surface facilitates suctioning and easy cleaning of the product

autoclavable, Material: PUR

Size	Color code	Length	Guedel Airway	Box
000	light blue	3.5 cm	REF 40-40-007	10
00	dark blue	5.0 cm	REF 40-40-008	10
0	black	6.0 cm	REF 40-40-009	10
1	white	7.0 cm	REF 40-40-010	10
2	green	8.5 cm	REF 40-40-020	10
3	yellow	9.0 cm	REF 40-40-030	10
4	red	10.0 cm	REF 40-40-040	10
5	purple	11.0 cm	REF 40-40-050	10



for single use, sterile
Material: EVA

Size	Color code	Length	Guedel Airway	Guedel Airway with elastic necktape	Box
000	light blue	3.5 cm	REF 40-20-007-1	REF 40-26-007-1	30
00	dark blue	5.0 cm	REF 40-20-008-1	REF 40-26-008-1	30
0	black	6.0 cm	REF 40-20-009-1	REF 40-26-009-1	30
1	white	7.0 cm	REF 40-20-010-1	REF 40-26-010-1	30
2	green	8.5 cm	REF 40-20-020-1	REF 40-26-020-1	30
3	yellow	9.0 cm	REF 40-20-030-1	REF 40-26-030-1	30
4	red	10.0 cm	REF 40-20-040-1	REF 40-26-040-1	30
5	purple	11.0 cm	REF 40-20-050-1	REF 40-26-050-1	30

BITEBLOCK

oral Tube fixations for ET Tube, Laryngeal Tube and Laryngeal Mask

- tube clips into biteblock and is therefore protected from occlusion



autoclavable, Material: PUR

Size	fits for	Biteblock	Box
4	ET Tube 6.5-8.0 mm	REF 40-70-040	10
5	ET Tube 8.5-11.0 mm	REF 40-70-050	10
	Laryngeal Tube # 2 + 2.5		
	Laryngeal Mask # 2 + 2.5		
6	Laryngeal Tube # 3, 4, 5	REF 40-70-060	10
	Laryngeal Mask # 3, 4, 5		

for single use, sterile
Material: EVA

Biteblock	Box	Biteblock with elastic necktape	Box
REF 40-60-040-1	25	REF 40-66-040-1	20
REF 40-60-050-1	25	REF 40-66-050-1	20
REF 40-60-060-1	25	REF 40-66-060-1	20