


Blood Exsanguination

Roll-On Cuffs

- Exsanguinator allows fast and easy blood exsanguination within seconds (as alternative to Esmarch Bandage)
- Material: Silicone, transparent
- Reusable, autoclavable up to 134°C (273°F)
- Latex free

Size	Color code	Patient limb	REF	
0	grey	16-21 cm	6"-8"	21-95-710
1	white	20-28 cm	8"-11"	21-95-711
2	light-blue	26-33 cm	10"-13"	21-95-712
3	yellow	31-40 cm	12"-16"	21-95-719
4	red	38-49 cm	15"-19"	21-95-722
5	green	47-60 cm	18"-23"	21-95-727
6	dark-blue	58-70 cm	23"-28"	21-95-729



Handpump
with manometer and
1 m (40") connecting tube
REF 21-95-350



Manometer
with 1 m (40") connecting tube
REF 21-95-300

Brake-Block

made of Silicone,
autoclavable up to 134°C (273°F)

small size **REF 21-95-200**
large size **REF 21-95-201**



Measuring Tape

color coded, for easy identification
of patient extremity and cuff size

REF 21-95-111



REF 21-95-400

Roll-On Cuff Set

small sizes

complete case with Roll-On Cuff # 1-4, manometer with
connecting tube, tape-measure, brake-block small and large

REF 21-95-400

Roll-On Cuff Set

large sizes

complete case with Roll-On Cuff # 3-6, handpump with
manometer and connecting tube, tape-measure, brake-block
small and large

REF 21-95-401

Esmarch Bandages

autoclavable up to 134°C (273°F), Box 1
 latex free
 Length: 5 m (16ft.)

Width	inch	REF
6 cm	2 ³ / ₈ "	20-20-060
8 cm	3 ¹ / ₈ "	20-20-080
10 cm	4"	20-20-100
12 cm	4 ³ / ₄ "	20-20-120
15 cm	6"	20-20-150



disposable, sterile packed, Box 10
 latex free

Length: 4 m (13 ft.)

Width	inch	REF
6 cm	2 ³ / ₈ "	20-21-060-1
8 cm	3 ¹ / ₈ "	20-21-080-1
10 cm	4"	20-21-100-1
15 cm	6"	20-21-150-1

Length: 3 m (9.8 ft.)

Width	inch	REF
8 cm	3 ¹ / ₈ "	20-23-080-1



Three Chamber Bloodexsanguination Garment

Automatic exsanguination for use with VBM Tourniquets 2800, 5000 and 5800

Due to the Three Chamber-System the garment is inflated automatically from distal to proximal.

With the 3 position zipper the garment can be adjusted to fit different size extremities for efficient exsanguination of the patients limb. Sleeve and bladder are exchangeable and available as spare parts.



REF 20-13-341



REF 20-13-343

Three Chamber Bloodexsanguination Garment, for arm
 Three Chamber Bloodexsanguination Garment, for leg

REF 20-13-341
 REF 20-13-343

