

Pressure Infusion Cuffs

500 ml and 1000 ml

Manual system for rapid pressure infusion and irrigation purposes.
Ideal for ambulances, operation rooms and intensive care stations.

500 ml

reusable

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

51-05-050

51-06-050

51-08-050



REF 51-08-050



REF 57-08-050

VBM Infusor

transparent Pressure Infusion Cuff
Wrap-Around Design, reusable

500 ml

with hand inflator and color coded pressure gauge

with hand inflator and fixed manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

57-02-050

57-06-050

57-08-050

1000 ml

pocket design with zip, reusable

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

52-05-100

52-06-100

52-08-100



REF 52-08-100

1000 ml

wrap around design, reusable

with hand inflator

with hand inflator and manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF

53-05-100

53-06-100

53-08-100



REF 53-08-100

for single patient use

500 ml

with hand inflator and pressure indicator

with hand inflator

REF

56-02-050

56-03-050

1000 ml

with hand inflator and pressure indicator

with hand inflator

REF

56-02-100

56-03-100



REF 56-02-100

VBM Infusor

transparent Pressure Infusion Cuff Wrap-Around Design, reusable

1000 ml

with hand inflator and color coded pressure gauge

with hand inflator and fixed manometer

with hand inflator and manometer
(shock resistant precision manometer)

REF


57-02-100

57-06-100

57-08-100



REF 57-08-100



Pressure Infusion Cuffs

1500 ml, 3000 ml and 5000 ml

For irrigation purposes e.g. during Arthroscopy, Laparoscopy and Hysteroscopy

1500 ml

pocket design with zip, reusable

with female Luer Lock

with hand inflator and manometer,
(shock resistant precision manometer)

with handpump and manometer
(shock resistant precision manometer)

REF

52-01-150

52-08-150

52-30-150



3000 ml

pocket design with zip, reusable

with female Luer Lock

with hand inflator and manometer,
(shock resistant precision manometer)

with handpump and manometer,
(shock resistant precision manometer)

REF

52-01-300

52-08-300

52-30-300

5000 ml

pocket design with zip, reusable

with female Luer Lock

with hand inflator and manometer,
(shock resistant precision manometer)

with handpump and manometer,
(shock resistant precision manometer)

REF

52-01-500

52-08-500

52-30-500



Accessories for Pressure Infusion Cuffs



Manometer 300 mmHg
shock resistant
REF 34-17-300



Manometer 300 mmHg
REF 20-16-300



Manometer 300 mmHg
with hand inflator and
release valve
REF 20-18-300



Manometer 300 mmHg
with hand inflator, release valve
and male Luer Lock
REF 51-09-666



Hand Inflator
with release valve
REF 51-10-666



Handpump
with release valve and
male Luer Lock
REF 51-30-000



Foot Pump
with 2 m (79") silicone hose,
release valve and male Luer Lock
REF 51-20-000



Hand Inflator
with pressure indicator, three-way
stop cock and male Luer Lock
REF 51-02-666



Hand Inflator
with three-way stop cock and
male Luer Lock
REF 51-03-666



Three-way stop cock
with two male and
one female Luer Lock
REF 51-02-000

