


# VBM Manual Suction Pump

An effective and powerful manual suction pump. For the quick suction of fluids from the oro- and nasopharyngeal cavities of adults and children in emergency situations.

- no power required
- portable and lightweight
- easy to squeeze with the maximum performance
- fits easily in most bags
- disposable collection canisters for adults and children
- cost effective alternative to electrical and foot powered devices
- soft and atraumatic suction tips



## Technical Data

Vacuum: max. 560mmHg  
size: 180 x 50 x 170mm  
weight: 170g

## Manual Suction Pump

complete with canister for adult and child, 3x suction tips (large, medium, small) and catheter adapter (for funnel shape suction catheters).

REF 30-55-500

# Accessories for Manual Suction Pump



**250ml Canister** for adults, with cap  
REF 30-56-250




**50ml Canister** for children, with cap  
REF 30-56-050

**Suction Tip, large, Ch 36**  
REF 30-57-812



15 O.D.

**Suction Tip, medium, Ch 28**  
REF 30-57-609



15 O.D.

**Suction Tip, small, Ch 14**  
REF 30-57-305



15 O.D.

15 I.D.



**Suction Cap** with 15mm I.D. connection  
REF 30-56-001

**Catheter Adapter** for funnel shape suction catheters  
REF 30-58-060

15 O.D.




**Carry Bag** fits pump, canisters, suction tips  
REF 30-64-000

**VBM** Medizintechnik GmbH

Einsteinstrasse 1 • DE-72172 Sulz a.N. • Germany

Tel.: +49 7454 / 95960 • Fax: +49 7454 / 959633 • E-mail: info@vbm-medical.de • www.vbm-medical.de



P325-1.0GB03.05