

VBM Endoscopy Mask

Allows simultaneous ventilation and endoscopy of the patient to increase patient comfort and the working quality of the doctor.

Application:

- **Fiberoptic Intubation**

The ET Tube can be railroaded orally or nasally along the bronchoscope through the flexible membrane. After intubation the membrane is pushed inside which allows easy removal of the mask.

- **Bronchoscopy**

The membrane guarantees an ideal seal and therefore an efficient ventilation before, during and after bronchoscopy without time pressure.

- **Gastroenterology**

- **Trans Esophageal Echocardiography**

The use of Entonox (50% N₂O, 50% O₂) is ideal for minor procedures where pain is predictable and of short duration. Entonox can be self administered from the conscious patient via the mask. It has a rapid analgesic effect and therefore allows a pain free insertion of the Endoscope through the hole in the membrane. Entonox wears off quickly without post operative problems.



Fiberoptic Intubation



ORDER INFORMATION

Endoscopy mask for Fiberoptic Intubation / Bronchoscopy

for single use, Material: PVC

Size	Patient	Membrane	REF	Box
0	Newborn	2.0 mm	30-40-000	6
1	Infant	2.0 mm	30-40-111	6
3	Child	3.0 mm	30-40-333	6
5	Adult	5.0 mm	30-40-555	6

Endoscopy mask for Gastroenterology / Trans Esophageal Echocardiography

for single use, Material: PVC

Size	Patient	Membrane	REF	Box
3	Child	5.0 mm	30-40-335	6
5	Adult	10.0 mm	30-40-777	6

FIXATION TAPES

Fixation tape

autoclavable, Material: Silicone

Patient	REF	Box
Child	35-70-155	10
Adult	35-70-255	10



Fixation tape

for single use, Material: Isoprene

Patient	REF	Box
Adult	35-30-255	10



Hook ring, Ø 32 mm
fits for Endoscopy masks,
Material: PA

REF 30-40-266	Box 10
----------------------	--------



Fixation tape

for single use, Material: Velours / Cotton

for Mask	REF	Box
# 1	35-40-622	10
# 3, 5	35-40-633	10



MEMBRANE

autoclavable, Material: Silicone

Hole	REF	Box
2.0 mm	30-40-121	5
3.0 mm	30-40-123	5
5.0 mm	30-40-125	5
10.0 mm	30-40-127	5



BRONCHOSCOPE AIRWAY

for single use, sterile, Material: EVA

The VBM Bronchoscope Airway avoids damage to the bronchoscope by the patient.

- the lumen is large enough to guide the bronchoscope and Endotracheal Tube
- the length allows to guide the bronchoscope into the pharynx
- the external biteblock protects the bronchoscope
- for single use - no risk of infection - no risk of contamination
- different sizes allow adjustment to patient and application

Size	Length	REF	Box
0	36 mm	30-40-400-1	10
2	65 mm	30-40-420-1	10
4	85 mm	30-40-440-1	10

