

VBM Stylets, Introducers & Tube Exchangers

METTS (Muallem ET Tube Stylet)

for difficult intubation

- core made of metal, malleable (maintains curvature)
- flexible pre-formed tip
- graduation marks
- for single use, sterile



METTI (Muallem ET Tube Introducer)

for difficult intubation and tube exchange

- core made of plastic, semi-rigid (flexible)
- flexible, pre-formed tip
- graduation marks
- for single use, sterile



INTRODUCER

for difficult intubation

- lumen-design, O₂-Connector included
- rigid pre-formed tip
- graduation marks
- for single use, sterile



ORDER INFORMATION

METTS

Size	Length	for ET Tube	REF	Box
8 CH	40 cm	≥ 3.5 mm	33-08-400-1	5
12 CH	65 cm	≥ 5.0 mm	33-12-650-1	5
14 CH	65 cm	≥ 6.0 mm	33-14-650-1	5

METTI

Size	Length	for ET Tube	REF	Box
12 CH	80 cm	≥ 5.0 mm	33-12-800-1	5
14 CH	80 cm	≥ 6.0 mm	33-14-800-1	5


Introducer

Size	Length	for ET Tube	REF	Box
15 CH	65 cm	≥ 6.0 mm	33-70-650-1	5

O₂-Connector

for oxygenation, fits all sizes with lumen-design (fits 11 CH, 14 CH, 15 CH, 19 CH), for single use

REF 33-60-005 Box 10



INTUBATION-CATHETER

for fiberoptically assisted intubation ($\leq 4,0$ mm)

- lumen-design, O₂-Connector included
- graduation marks
- for single use, sterile



TUBE EXCHANGER

for extubation and exchange of Endotracheal Tubes

- lumen-design, O₂-Connector included
- graduation marks
- for single use, sterile



STYLETS

malleable intubation stylet

- soft, atraumatic tip
- with and without moveable silicone connector
- autoclavable



ORDER INFORMATION

Intubation-Catheter

Size	Length	for ET Tube	for Fiberoptik	REF	Box
19 CH	56 cm	≥ 7.0 mm	≤ 4.0 mm	33-60-560-1	5

Tube Exchanger

Size	Length	for ET Tube	REF	Box
11 CH	80 cm	≥ 4.0 mm	33-45-800-1	5
14 CH	80 cm	≥ 5.0 mm	33-46-800-1	5
19 CH	80 cm	≥ 7.0 mm	33-47-800-1	5

Stylet

Size	Length	for ET Tube	without silicone connector	Box	with moveable silicone connector (15 mm I.D.)	Box
					REF	
small	34 cm	2.5 - 3.5 mm	REF 90-10-111	25	REF 90-10-001	10
medium	38 cm	4.0 - 4.5 mm	REF 90-10-222	25	REF 90-10-002	10
large	45 cm	5.0 - 6.5 mm	REF 90-10-333	25	REF 90-10-003	10
large	45 cm	7.0 - 11.0 mm	REF 90-10-444	25	REF 90-10-004	10